


Product Datasheet

Rat Fgfbp1 ELISA Kit (orb1292631)

Description	Rat Fgfbp1 (Fibroblast growth factor-binding protein 1) ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	Fgfbp1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9QY10
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat Fgfbp1(Fibroblast growth factor-binding...