

Product Datasheet

Rat Htr1a ELISA Kit (orb1292643)

Description Rat Htr1a (5-hydroxytryptamine receptor 1A)ELISA Kit

Reactivity Rat

Range 78.125-5000pg/ml

Target Htr1a

Storage 2-8 °C for 6 months


Note For research use only

Sample Types Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples

Uniprot ID **P19327**

Sensitivity 46.875pg/ml

Expiration Date Please enquire.



Rat Htr1a(5-hydroxytryptamine receptor 1...