


Product Datasheet

Rat 8-OHdG ELISA Kit (orb1292675)

Description	Rat 8-OHdG (8-Hydroxydeoxyguanosine)
Reactivity	Rat
Range	0.469-30ng/ml
Target	8-OHdG
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, urine, Other liquid samples
Sensitivity	0.281ng/ml
Expiration Date	Please enquire.



Rat 8-OHdG(8-Hydroxydeoxyguanosine)