


## Product Datasheet

### Porcine HS ELISA Kit (orb1292682)

<b>Description</b>	Porcine HS (Heparan Sulfate) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	31.25-2000ng/ml
<b>Target</b>	HS
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Sensitivity</b>	18.75ng/ml
<b>Expiration Date</b>	Please enquire.



Porcine  
HS(Heparan  
Sulfate)  
ELISA Kit