

Product Datasheet

Porcine GPX1 ELISA Kit (orb1292694)

Description Porcine GPX1 (Glutathione peroxidase 1) ELISA Kit

Reactivity Porcine

Range 0.156-10ng/ml

Target GPX1

Storage 2-8 °C for 6 months


Note For research use only

Sample Types Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples

Uniprot ID **Q8MJ14**

Sensitivity 0.094ng/ml

Expiration Date Please enquire.



Porcine
GPX1(Glutathione
peroxidase 1) E...