


## Product Datasheet

### Monkey INS ELISA Kit (orb1292708)

<b>Description</b>	Monkey INS (Insulin) ELISA Kit
<b>Reactivity</b>	Monkey
<b>Range</b>	7.813-500pg/ml [1ng=10uIU]
<b>Target</b>	INS
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>P30407</b>
<b>Sensitivity</b>	4.688pg/ml
<b>Expiration Date</b>	Please enquire.



Monkey  
INS(Insulin)  
ELISA Kit