


## Product Datasheet

### Mouse Tyrobp ELISA Kit (orb1292723)

<b>Description</b>	Mouse Tyrobp (TYRO protein tyrosine kinase-binding protein) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	Tyrobp
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>O54885</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
Tyrobp  
(TYRO  
protein  
tyrosine  
kina...