


## Product Datasheet

### Mouse Mboat4 ELISA Kit (orb1292732)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Mouse Mboat4 ELISA Kit   |
| <b>Reactivity</b>      | Mouse  |
| <b>Range</b>           | 78.125-5000pg/ml   |
| <b>Target</b>          | Mboat4   |
| <b>Storage</b>         | 2-8 °C for 6 months  |
| <b>Note</b>            | For research use only  |
| <b>Sample Types</b>    | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| <b>Uniprot ID</b>      | <b>POC7A3</b>  |
| <b>Sensitivity</b>     | 46.875pg/ml  |
| <b>Expiration Date</b> | Please enquire.  |



Mouse Mboat4  
(Ghrelin O-  
acyltransferase)...