


Product Datasheet

Mouse Icos ELISA Kit (orb1292738)

Description	Mouse Icos ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/ml
Target	Icos
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9WVS0
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse Icos
(Inducible T-
cell
costimulato...