


Product Datasheet

Mouse Hsd17b3 ELISA Kit (orb1292739)

Description	Mouse Hsd17b3 ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	Hsd17b3
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P70385
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse
Hsd17b3
(Testosterone
17-beta-
dehy...