


Product Datasheet

Mouse Ddr1 ELISA Kit (orb1292747)

Description	Mouse Ddr1 ELISA Kit
Reactivity	Mouse
Range	15.62-1000pg/ml
Target	Ddr1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q03146
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse Ddr1
(Epithelial discoidin domain-...