


Product Datasheet

Mouse Pomgnt1 ELISA Kit (orb1292762)

Description	Mouse Pomgnt1 (Protein O-linked-mannose beta-1,2-N-acetylglucosaminyltransferase 1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Pomgnt1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Cell or tissue lysate, Other liquid samples
Uniprot ID	Q91X88
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
Pomgnt1(Protein
O-linked-
mannose b...