


Product Datasheet

Mouse Prdx6 ELISA Kit (orb1292765)

| | |
|------------------------|--|
| Description | Mouse Prdx6 (Peroxiredoxin-6) ELISA Kit |
| Reactivity | Mouse |
| Range | 15.625-1000pg/ml |
| Target | Prdx6 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Sample Types | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| Uniprot ID | 008709 |
| Sensitivity | 9.375pg/ml |
| Expiration Date | Please enquire. |



Mouse
Prdx6(Peroxiredoxin-
6) ELISA Kit