


Product Datasheet

Mouse Efna1 ELISA Kit (orb1292782)

Description	Mouse Efna1 (Ephrin-A1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Efna1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P52793
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
Efna1(Ephrin-
A1) ELISA Kit