


Product Datasheet

Mouse Ppbp ELISA Kit (orb1292787)

Description	Mouse Ppbp/Cxcl7 (Chemokine subfamily B Cys-X-Cys) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	Ppbp/Cxcl7
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Gene ID:57349
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
Ppbp/Cxcl7(Chemokine
subfamily B C...