


Product Datasheet

Mouse Wfdc12 ELISA Kit (orb1292789)

Description	Mouse Wfdc12 ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	Wfdc12
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9JHY3
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
Wfdc12
(Single WAP
motif protein
2...