


Product Datasheet

Mouse Sstr2 ELISA Kit (orb1292813)

Description	Mouse Sstr2 (Somatostatin receptor type 2)ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	Sstr2
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P30875
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse Sstr2
(Somatostatin
receptor type
...)