


## Product Datasheet

### Mouse Sfrp4 ELISA Kit (orb1292854)

<b>Description</b>	Mouse Sfrp4 (Secreted frizzled-related protein 4) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	Sfrp4
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>Q9Z1N6</b>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse Sfrp4(Secreted frizzled-related pr...