


## Product Datasheet

### Human MBD2 ELISA Kit (orb1292864)

<b>Description</b>	Human MBD2 (Methyl-CpG-binding domain protein 2) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	MBD2
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
MBD2(Methyl-  
CpG-binding  
domain pro...