


## Product Datasheet

### Human HTR3A ELISA Kit (orb1292896)

<b>Description</b>	Human HTR3A (5-hydroxytryptamine receptor 3A) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	HTR3A
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P46098</b>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Human HTR3A(5-hydroxytryptamine receptor...