


Product Datasheet

Human GALNT14 ELISA Kit (orb1292926)

Description	Human GALNT14 (Polypeptide N-acetylgalactosaminyltransferase 14) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	GALNT14
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q96FL9
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
GALNT14(Polypeptide
N-acetylgalact...