


Product Datasheet

Human CHRNA3 ELISA Kit (orb1292969)

Description	Human CHRNA3 (Neuronal acetylcholine receptor subunit alpha-3) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	CHRNA3
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P32297
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
CHRNA3(Neuronal
acetylcholine
rece...