


Product Datasheet

Human BNIP3L ELISA Kit (orb1292981)

Description	Human BNIP3L (BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	BNIP3L
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O60238
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
BNIP3L(BCL2/adenovirus
E1B 19 kDa ...