


## Product Datasheet

### Human UBE2R2 ELISA Kit (orb1293052)

<b>Description</b>	Human UBE2R2 (Ubiquitin-conjugating enzyme E2 R2) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	UBE2R2
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q712K3</b>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
UBE2R2(Ubiquitin-  
conjugating  
enzym...