


## Product Datasheet

### Human TAU-1 ELISA Kit (orb1293087)

<b>Description</b>	Human TAU-1 (microtubule-associated protein tau-1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	TAU
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P10636</b>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Human TAU-1 (microtubule-associated prote...