


Product Datasheet

Human ADAM 15 ELISA Kit (orb1293091)

Description	Human ADAM 15 (Disintegrin and metalloproteinase domain-containing protein 15) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	ADAM 15
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q13444
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human ADAM 15(Disintegrin and metallopro...