


## Product Datasheet

### Human WDR5 ELISA Kit (orb1293093)

<b>Description</b>	Human WDR5 (WD repeat domain 5)ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	WDR5
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P61964</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
WDR5 (WD  
repeat  
domain 5)  
ELISA Ki...