


Product Datasheet

Human PGAM5 ELISA Kit (orb1293130)

| | |
|------------------------|---|
| Description | Human PGAM5 (Serine/threonine-protein phosphatase PGAM5, mitochondrial) ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/ml |
| Target | PGAM5 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | Q96HS1 |
| Sensitivity | 0.094ng/ml |
| Expiration Date | Please enquire. |



Human PGAM5
(Serine/threonine-
protein ph...