


## Product Datasheet

### Human GPR37L1 ELISA Kit (orb1293159)

<b>Description</b>	Human GPR37L1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	GPR37L1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>O60883</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
GPR37L1  
(Prosaposin  
receptor  
GPR37...