


## Product Datasheet

### Human SNAI1 ELISA Kit (orb1293239)

<b>Description</b>	Human SNAI1 (Snail family transcriptional repressor 1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	SNAI1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>O95863</b>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
SNAI1(Snail  
family  
transcriptional...