


Product Datasheet

Human SIGLEC5 ELISA Kit (orb1293245)

Description	Human SIGLEC5 (Sialic acid-binding Ig-like lectin 5) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	SIGLEC5
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O15389
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
SIGLEC5(Sialic
acid-binding
Ig-lik...