


Product Datasheet

Human PVRL4 ELISA Kit (orb1293269)

Description	Human PVRL4/Nectin-4 (Nectin cell adhesion molecule 4) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	Nectin-4
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q96NY8
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
PVRL4/Nectin-
4(Nectin cell
adhesio...