


## Product Datasheet

### Human PKNOX2 ELISA Kit (orb1293285)

<b>Description</b>	Human PKNOX2 (Homeobox protein PKNOX2) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78.125-5000pg/ml
<b>Target</b>	PKNOX2
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q96KN3</b>
<b>Sensitivity</b>	46.875pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
PKNOX2(Homeobox  
protein PKNOX2)  
EL...