


## Product Datasheet

### Human PCYOX1 ELISA Kit (orb1293296)

<b>Description</b>	Human PCYOX1 (Prenylcysteine oxidase 1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	PCYOX1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9UHG3</b>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
PCYOX1(Prenylcysteine  
oxidase 1) E...