


## Product Datasheet

### Guinea pig IL-18 ELISA Kit (orb1293322)

<b>Description</b>	Guinea pig IL-18 (Interleukin 18) ELISA Kit
<b>Reactivity</b>	Guinea pig
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	IL-18
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Guinea pig IL-18(Interleukin 18) ELISA K...