


Product Datasheet

Bovine IL1ra ELISA Kit (orb1293348)

Description	Bovine IL1ra (Interleukin 1 Receptor Antagonist) ELISA Kit
Reactivity	Bovine
Range	31.25-2000pg/ml
Target	IL-1ra
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	077482
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Bovine
IL1ra(Interleukin
1 Receptor
Anta...