



Product Datasheet

Occludin antibody (orb1294255)


Description	Rabbit polyclonal antibody to Occludin
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / Extracellular domain corresponding to Human Occludin / OCLN
Target	Occludin
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 1% BSA, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500-1:1500 Immunofluorescence 1:200-1:300 Immunohistochemistry 1:100
Expiration Date	12 months from date of receipt.



Anti-Occludin antibody - at 1/1000 dilut...



Anti-Occludin antibody - at 1/1000 dilut...



Immunohistochemical analysis of paraffin...