


Product Datasheet

Enterovirus 71 antibody (orb1294272)

Description	Rabbit polyclonal antibody to Enterovirus 71
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fragment corresponding to Enterovirus 71 / EV71
Target	Enterovirus 71 / EV71
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:5000-1:10000 Immunofluorescence 1:1000 - 1:2000
Expiration Date	12 months from date of receipt.



Enterovirus 71 / EV71 antibody at 1/5000 ...



Specificity and sensitivity of Anti Ente...