


## Product Datasheet

### GST tag [6D] monoclonal Antibody (orb1294274)

<b>Description</b>	Mouse Monoclonal antibody to GST-tag (6D) Monoclonal
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Synthetic peptide
<b>Target</b>	GST tag [6D] monoclonal
<b>Storage</b>	Store at +4°C for short term storage. Long time storage is recommended at -20°C 50mM Tris Glycine, 40% Glycerol (pH7.4). 0.025% ProClin 300 was added as a preservative.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	Western Blot 1:2000-1:5000 ELISA 1:5000-1:10000
<b>Expiration Date</b>	12 months from date of receipt.



All lanes:  
Anti GST-tag [6D] monoclonal

...