


Product Datasheet

His tag Antibody (orb1294276)

Description	Rabbit polyclonal antibody to His tag
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	KLH-conjugated HHHHHH (6X His)
Target	His tag
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:2000-1:4000 Immunofluorescence 1:500-1:1000
Expiration Date	12 months from date of receipt.



His-tag antibody detects His-tagged prot...