


Product Datasheet

HA tag Antibody (orb1294277)

Description	Rabbit polyclonal antibody to HA tag
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Synthetic peptide (YPYDVPDYA)
Target	HA tag
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:3000-1:7000 Immunofluorescence 1:500-1:1000
Expiration Date	12 months from date of receipt.



All lanes: Anti HA-tag antibody at 1/500...