



Product Datasheet

SARS-CoV-2 (COVID-19) Spike Antibody (orb1294278)

Description	Rabbit polyclonal antibody to SARS CoV 2 spike glycoprotein
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IF, NeA, WB
Immunogen	Recombinant fragment corresponding to SARS CoV 2 spike glycoprotein
Target	SARS-CoV-2 (COVID-19) Spike
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:5000-1:18000 Immunofluorescence 1:500-1:1000
Expiration Date	12 months from date of receipt.



Anti-SARS-CoV-2 spike protein antibody a...



Immunofluorescence: Cells were stained W...