


Product Datasheet


SARS-CoV-2 (COVID-19) nucleocapsid Antibody (orb1294283)

Description	Rabbit polyclonal antibody to SARS-CoV 2 nucleocapsid protein
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IF, NeA, WB
Immunogen	Recombinant fragment corresponding to SARS-CoV 2 nucleocapsid protein
Target	SARS-CoV-2 (COVID-19) nucleocapsid
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:3000-1:5000 Immunofluorescence 1:500-1:1000
Expiration Date	12 months from date of receipt.



All lanes: Anti-SARS-CoV-2 nucleocapsid

...



Immunofluorescence:

Cells were stained

W...