



Product Datasheet

CTGF Antibody (orb1294294)


Description	Rabbit polyclonal antibody to CTGF
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / encompassing a sequence within the center region
Target	CTGF
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 40% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500-1:1000 Immunofluorescence 1:200-1:300 Immunohistochemistry (Paraffin) 1: 100-1:150
Expiration Date	12 months from date of receipt.



Anti-CTGF antibody at 1/700 dilution, Ly...



Anti-CTGF antibody at 1/700 dilution, Ly...



Immunofluorescence: cells were fixed wit...