



Product Datasheet

Survivin antibody (orb1294349)


Description	Rabbit polyclonal antibody to Survivin
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / encompassing a sequence within the C-terminus region.
Target	Survivin / BIRC5
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:2000-1:3000 Immunofluorescence 1:400-1:800 Immunohistochemistry (Paraffin) 1:100-1:300
Expiration Date	12 months from date of receipt.



Survivin/BIRC5 antibody at 1/2000 diluti...



Cells were fixed with 4% PFA for 10 min ...



Immunohistochemical analysis of paraffin...