



Product Datasheet

E-cadherin (Extracellular domain) antibody (orb1294382)


Description	Rabbit polyclonal antibody to E-cadherin antibody (Extracellular domain)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / Extracellular domain corresponding to Human E-cadherin
Target	E-cadherin (Extracellular)
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500-1:1000 Immunofluorescence 1:200-1:500 Immunohistochemistry (Paraffin) 1:100-1:200
Expiration Date	12 months from date of receipt.



Anti-E cadherin antibody at 1/500 diluti...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...