



## Product Datasheet

# NOTCH1 Antibody (orb1294410)


<b>Description</b>	Rabbit polyclonal antibody to NOTCH1
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, WB
<b>Immunogen</b>	Synthetic peptide / encompassing a sequence within the center region.
<b>Target</b>	NOTCH1
<b>Storage</b>	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 1% BSA, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit-IgG
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	Western Blot 1:500 Immunofluorescence 1:200 - 1:500
<b>Expiration Date</b>	12 months from date of receipt.



All lanes: Anti-  
NOTCH1  
antibody at  
1/700...



All lanes: Anti-  
NOTCH1  
antibody at  
1/700...



Sample:  
Neuron 2A;  
NOTCH1  
(orb1294410):

...