


Product Datasheet

EIF4EBP1 (phospho-Ser65) antibody (orb129595)

Description	Rabbit polyclonal to EIF4EBP1 (phospho-Ser65).
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphorylated peptide around S65 of human Phospho-eIF4EBP1-S65(NP_004086.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20KDa
Uniprot ID	Q13541
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis on 293 cell lysate...