


## Product Datasheet

### Flag Tag antibody (orb129610)

<b>Description</b>	Mouse monoclonal to Flag Tag.FLAG tag Mouse mAb is part of the series of Tag antibodies, the excelle...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Other
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IF, IHC-P, IP, WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	WB: 1:3000-1:10000, IF:1:200
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK 293 cell ly...