


Product Datasheet

Phospho-INSR-Y1150/IGF1R-Y1135 pAb (orb129649)

Description	Rabbit polyclonal antibody to Phospho-INSR-Y1150/IGF1R-Y1135 pAb
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding Y1150 of human INSR and Y1135 of human IGF1R
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	80kDa
Uniprot ID	P08069
Entrez	34803643
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...